

VBM Pressure Infusion Cuff

- transparent -

Transparent

- allows clear view on the fluid bag and the fluid level

Shock Resistant Precision Manometer

- exact pressure control
- overpressure proved
- shock resistant and long-lasting

Wrap-Around Design

- fluid bag can be replaced quickly and easily
- consistent pressure on the complete cuff

PU Material

- easy to clean
- reusable
- durable
- latex free




Stop Cock

- allows quick and complete deflation of the cuff

Large Hand Inflator

- for quick inflation of the cuff



**Pressure Infusion Cuff 500ml,
with Manometer and Hand Inflator**

REF 57-08-050

**Pressure Infusion Cuff 1000ml,
with Manometer and Hand Inflator**

REF 57-08-100



**Pressure Infusion Cuff 500ml,
with Hand Inflator and Pressure Indicator**

REF 57-02-050

**Pressure Infusion Cuff 1000ml,
with Hand Inflator and Pressure Indicator**

REF 57-02-100



Components



Cuff 500ml

REF 57-71-050

Cuff 1000ml

REF 57-71-100



Manometer 0-300mmHg

REF 34-17-300



Hand Inflator

with Pressure Indicator

REF 51-02-666



Hand Inflator

REF 51-03-666

VBM Medizintechnik GmbH

Einsteinstrasse 1 • DE-72172 Sulz a.N. • Germany

Tel.: +49 7454 / 959610 • Fax: +49 7454 / 959633 • e-mail: info@vbm-medical.de • www.vbm-medical.de

