

Pressure Infusion Cuffs

Pressure Infusion Cuffs

Manual system for rapid I.V. infusion of blood and solutions. Ideal for ambulances, operating rooms and intensive care stations



Pressure Infusion Cuffs 500 ml reusable

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

51-05-050

51-06-050

51-08-050



Pressure Infusion Cuffs 1000 ml pocket version with zip, reusable

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

52-05-100

52-06-100

52-08-100



Pressure Infusion Cuffs 1000 ml Wrap-Around Design with velcro, reusable

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

53-05-100

53-06-100

53-08-100



for single patient use

500 ml

with hand inflator, three-way stop cock and pressure indicator
with hand inflator and three-way stop cock

REF

56-02-050

56-03-050

1000 ml

with hand inflator, three-way stop cock and pressure indicator
with hand inflator and three-way stop cock

REF

56-02-100

56-03-100



VBM Infusor

Transparent Pressure Infusion Cuffs
Wrap-Around Design, reusable

500 ml

with hand inflator and color coded pressure indicator

REF

57-02-050

with manometer and hand inflator
(shock resistant precision manometer)

57-08-050

1000 ml

with hand inflator and color coded pressure indicator

REF

57-02-100

with manometer and hand inflator
(shock resistant precision manometer)

57-08-100

Tourniquet "Tourny I"

with stainless steel buckle

REF 20-30-000 Box 20



Tourniquet "Tourny II"

with aluminum buckle

REF 20-32-000 Box 20



Tourniquet "Tourny III"

with plastic buckle

REF 20-30-100 Box 10



Tourniquet for Child

elastic strap with velcro fastening


REF 20-30-005 Box 10



Tourniquet Strap

elastic strap 40 cm,
25 pcs perforated on one roll

REF 20-30-400 Box 1 roll



Manual Suction Pump

An effective and powerful manual suction pump for the quick suction of fluids from the oro- and nasopharyngeal cavities in emergency situations.

Complete Set

with canister for adult and child, 3x suction tips (large, medium, small) and catheter adapter (for funnel shape suction catheters). color: blue

REF 30-55-500

ditto but, color: olive (military)

REF 30-55-501

- portable and lightweight
- always ready for use as independent from any energy source
- easy to squeeze with maximum performance
- fits easily in most bags
- disposable collection canisters for adults and children
- cost effective alternative to electrical and foot powered devices
- soft and atraumatic suction tips
- valve prevents backflow of liquids to the patient
- MRI proof
- 15 mm I.D. connection
 - fits directly to the ET Tube
 - allows connection of included suction tips
 - allows connection of included catheter adapter for funnel shape suction catheters



250 ml Canister

for adults, with cap

REF 30-56-250 Box 10



50 ml Canister

for children, with cap

REF 30-56-050 Box 10

Suction Cap

with 15 mm I.D. connection

REF 30-56-001 Box 10



Catheter Adapter

for funnel shape suction catheters

REF 30-58-060 Box 10



Carry Bag

fits pump, canisters, suction tips

REF 30-64-000



Suction Tip

large, 36 CH

REF 30-57-812 Box 10



medium, 28 CH

REF 30-57-609 Box 10



small, 14 CH

REF 30-57-305 Box 10

