


## Guedel Airways

### oropharyngeal airway with ET Tube fixation possibility

- rail and lateral clamp design allow direct fixation of the ET Tube to the Guedel Airway
- hard biteblock but soft and atraumatic tube
- same lumen from proximal to distal end allows the use of largest possible suction catheters
- one side open with smooth inner surface facilitates suctioning and easy cleaning of the product



Size	Color code	Length (cm)	autoclavable up to 134°C (273°F) Material: PU	Box	for single use, sterile Material: EVA	Box	for single use, sterile, with elastic necktape Material: EVA	Box
000	light blue	3.5	40-40-007	10	40-20-007-1	30	40-26-007-1	30
00	marine	5.0	40-40-008	10	40-20-008-1	30	40-26-008-1	30
0	black	6.0	40-40-009	10	40-20-009-1	30	40-26-009-1	30
1	white	7.0	40-40-010	10	40-20-010-1	30	40-26-010-1	30
2	green	8.5	40-40-020	10	40-20-020-1	30	40-26-020-1	30
3	yellow	9.0	40-40-030	10	40-20-030-1	30	40-26-030-1	30
4	red	10.0	40-40-040	10	40-20-040-1	30	40-26-040-1	30
5	purple	11.0	40-40-050	10	40-20-050-1	30	40-26-050-1	30

## Biteblocks

### oral Tube fixations for ET Tube, Laryngeal Tube and Laryngeal Mask

- tube clips into biteblock and is therefore protected from occlusion



Size	fits for	Size	autoclavable up to 134°C Material: PU	Box	for single use, sterile Material: EVA	Box	for single use, sterile, with elastic necktape Material: EVA	Box
4	ET Tube	6.5-8.0 mm	40-70-040	10	40-60-040-1	25	40-66-040-1	20
5	ET Tube	8.5-11.0 mm	40-70-050	10	40-60-050-1	25	40-66-050-1	20
	Laryngeal Tube	2 + 2.5						
	Laryngeal Mask	2 + 2.5						
6	Laryngeal Tube	3, 4, 5	40-70-060	10	40-60-060-1	25	40-66-060-1	20
	Laryngeal Mask	3, 4, 5						

## Fixation Tapes for Guedel Airways and Biteblocks

### Foam Necktape

REF 40-25-000

Box 5



### Elastic Necktapes

Length	REF	Box
40 cm	40-06-040	10
55 cm	40-06-055	10
75 cm	40-06-075	10



### Foam Pads

self-adhesive, for Airway and Biteblock

Size	For size	REF	Box
small	000, 1, 2	40-20-011	30
large	00, 0, 3, 4, 5	40-20-033	30

