

Masks

Endoscopy Masks

Allow simultaneous ventilation and endoscopy of the patient to increase patient comfort and the working quality of the doctor.

for single use

Application:

- Fiberoptic Intubation
The ET tube can be railroaded orally or nasally along the bronchoscope through the flexible membrane. After intubation the membrane is pushed inside which allows easy removal of the mask.
- Bronchoscopy
The membrane guarantees an ideal seal and therefore an efficient ventilation before, during and after bronchoscopy without time pressure.
- Gastroenterology
- Trans Esophageal Echocardiography

Size	Patient	Membrane	REF	Box
0	Newborn	2.0 mm	30-40-000	6
1	Infant	2.0 mm	30-40-111	6
3	Child	3.0 mm	30-40-333	6
3	Child	5.0 mm	30-40-335	6
5	Adult	5.0 mm	30-40-555	6
5	Adult	10.0 mm	30-40-777	6



Harness

Material: Silicone, latex free reusable, autoclavable up to 134°C

Patient	REF	Box
Child	35-70-155	10
Adult	35-70-255	10



Fixation Tapes

Material: foam

For Mask	REF	Box
# 1	35-40-622	10
# 3+5	35-40-633	10



Hook Ring

made of PE, Ø 32 mm fits for Endoscopy masks

REF 30-40-266 Box 10



Bronchoscope Airways

with external biteblock, single use, sterile

- protect expensive equipment from damage
- guide bronchoscope and ET Tube

Size	Length	REF	Box
0	36 mm	30-40-400-1	10
2	65 mm	30-40-420-1	10
4	85 mm	30-40-440-1	10



Face Masks made of Silicone

latex-free, reusable, autoclavable up to 134°C

with plastic dome and inflatable cushion

Size	Patient	Color	REF	Box
2	Child	green	35-70-222	1
3	Adult small	yellow	35-70-223	1
4	Adult medium	red	35-70-224	1
5	Adult large	purple	35-70-225	1
6	Adult extra large	blue	35-70-226	1



with plastic dome

Size	Patient	Color	REF	Box
3	Adult small	yellow	80-11-003	10
4	Adult medium	red	80-11-004	10
5	Adult large	purple	80-11-005	10



Bi-Masks

Silicone mask with internal plastic shell (PSU)

Size	Patient	Color	REF	Box
3	Adult small	yellow	35-65-223	1
4	Adult medium	red	35-65-224	1
5	Adult large	purple	35-65-225	1
6	Adult extra large	blue	35-65-226	1



Masks for Children, completely made of Silicone

latex-free, autoclavable up to 134°C

Rendell Baker Masks

Size	Patient	REF	Box
0	Neonate	35-60-000	10
1	Infant	35-60-001	10
2	Child small	35-60-002	10
3	Child	35-60-003	10



Masks, round

Size	Patient	REF	Box
0	Neonate	80-11-000	10
1	Infant	80-11-001	10
2	Child	80-11-002	10



Resuscitation Bags, single use

PVC Face Masks

with inflatable cushion, for single use


Size	Patient	Color code	Connection	REF	Box
1	Neonate	-	15 mm O.D.	35-40-201	10
2	Infant	-	15 mm O.D.	35-40-202	10
3	Child	yellow	22 mm I.D.	35-40-203	10
4	Adult small	red	22 mm I.D.	35-40-204	10
5	Adult medium	purple	22 mm I.D.	35-40-205	10
6	Adult large	blue	22 mm I.D.	35-40-206	10



Harness

latex free, for single use

Patient	Material	REF	Box
Adult	synthetic rubber	35-30-255	10



Fixation Tapes

for face masks for single use
material: foam

For mask	Length	REF	Box
# 1	350 mm	35-40-610	10
# 2-3	500 mm	35-40-611	10
# 4-6	600 mm	35-40-633	10



PVC Resuscitation Bags

for single use

Adult 1800 ml	REF	Box
with 40 cmH ₂ O pressure relief	85-10-100	1
without pressure relief	85-10-101	1



Child 500 ml	REF	Box
with 40 cmH ₂ O pressure relief	85-10-200	1
without pressure relief	85-10-201	1



Infant 250 ml	REF	Box
with 40 cmH ₂ O pressure relief	85-10-300	1
without pressure relief	85-10-301	1



Masks

CPAP Respiration Sets

single patient use

consisting of:

CPAP Mask, Harness, PEEP Valve, adjustable from 5 - 20 cmH₂O and Adapter for 22 mm tube connection

Size	REF	Box
medium	35-45-001	1
large	35-45-002	1

ditto, but with PEEP Valve, adjustable from 2.5 - 10 cmH₂O

Size	REF	Box
medium	35-45-003	1
large	35-45-004	1



CPAP Masks

with inflatable cushion and with Adapter REF 60-50-007

Size	REF	Box
medium	35-40-100	5
large	35-40-000	5



PEEP Valves

for single use, MRI-proof

adjustable from	REF	Box
5 - 20 cmH ₂ O, blue	85-13-002	1
2.5 - 10 cmH ₂ O, green	85-13-004	1



30 I.D.

PEEP Valves

autoclavable up to 134°C, MRI-proof

adjustable from	REF	Box
5 - 20 cmH ₂ O, blue	80-13-002	1
2.5 - 10 cmH ₂ O, green	80-13-004	1




30 I.D.

Harness

latex free, for single use

REF 35-30-255 Box 10



Adapter

for single use

REF 60-50-007 Box 10



29 I.D./
32 O.D.

22 O.D.