

Tube Inhaler


for application of MDI canisters on intubated patients

autoclavable up to 134°C, Material: PSU

- compatible with any MDI canister
- spray application during inspiration
- easy to use
- no interruption of ventilation
- suction port (max. 16 CH)




Suctioning with max. 16 CH



Tube Inhaler


Catheter Mount Tubings made of Silicone with Tube Inhaler
 autoclavable up to 134°C

with Tube Inhaler „Single Swivel“ with suction port



Length: 13 cm	16 cm	20 cm	23 cm	Box
65-10-009		65-10-029		10
65-30-009		65-30-029		10
	65-40-009		64-40-029	10

with Tube Inhaler „Double Swivel“ with suction port



Length: 13 cm	16 cm	20 cm	23 cm	Box
65-15-209		65-15-229		10
65-35-209		65-35-229		10
	65-45-209		64-45-229	10

Suction cap

made of Silicone, clear,
 suitable for Tube Inhaler

REF 60-60-119 Box 10

