


Tubings - smooth / clear

reusable, autoclavable up to 134°C (273°F)

Material: Silicone (SI)

Box 1 roll = 25 m (82ft.)

I.D. (mm)	x	C	O.D. (mm)	I.D. (inch)	x	C	O.D. (inch)	REF
1.0	x	0.5	2.0	.04"	x	.02"	.08"	68-010-05
1.0	x	1.0	3.0	.04"	x	.05"	.12"	68-010-10
1.5	x	0.75	3.0	.06"	x	.03"	.12"	68-015-075
2.0	x	1.0	4.0	.08"	x	.05"	.16"	68-020-10
3.0	x	1.0	5.0	.12"	x	.05"	.20"	68-030-10
3.0	x	1.5	6.0	.12"	x	.06"	.24"	68-030-15
4.0	x	1.0	6.0	.16"	x	.05"	.24"	68-040-10
4.0	x	1.5	7.0	.16"	x	.06"	.28"	68-040-15
4.0	x	2.0	8.0	.16"	x	.08"	.31"	68-040-20
5.0	x	1.0	7.0	.20"	x	.05"	.28"	68-050-10
5.0	x	1.5	8.0	.20"	x	.06"	.31"	68-050-15
5.0	x	2.0	9.0	.20"	x	.08"	.35"	68-050-20
6.0	x	1.0	8.0	.24"	x	.05"	.31"	68-060-10
6.0	x	1.5	9.0	.24"	x	.06"	.35"	68-060-15
6.0	x	2.0	10.0	.24"	x	.08"	.39"	68-060-20
6.0	x	3.0	12.0	.24"	x	.12"	.47"	68-060-30
7.0	x	1.5	10.0	.28"	x	.06"	.39"	68-070-15
7.0	x	2.0	11.0	.28"	x	.08"	.43"	68-070-20
7.0	x	3.0	13.0	.28"	x	.12"	.51"	68-070-30
8.0	x	1.0	10.0	.31"	x	.05"	.39"	68-080-10
8.0	x	2.0	12.0	.31"	x	.08"	.47"	68-080-20
8.0	x	3.0	14.0	.31"	x	.12"	.55"	68-080-30
9.0	x	2.0	13.0	.35"	x	.08"	.51"	68-090-20
10.0	x	2.0	14.0	.39"	x	.08"	.55"	68-100-20
10.0	x	3.0	16.0	.39"	x	.12"	.63"	68-100-30
12.0	x	3.0	18.0	.47"	x	.12"	.71"	68-120-30
12.0	x	4.0	20.0	.47"	x	.16"	.79"	68-120-40
14.0	x	3.0	20.0	.55"	x	.12"	.79"	68-140-30
15.0	x	3.0	21.0	.59"	x	.12"	.83"	68-150-30
16.0	x	3.0	22.0	.63"	x	.12"	.87"	68-160-30
18.0	x	3.0	24.0	.71"	x	.12"	.94"	68-180-30
20.0	x	2.0	24.0	.79"	x	.08"	.94"	68-200-20
20.0	x	3.0	26.0	.79"	x	.12"	1.02"	68-200-30
25.0	x	2.0	29.0	.98"	x	.08"	1.14"	68-250-20
25.0	x	4.0	33.0	.98"	x	.16"	1.30"	68-250-40



Further sizes are available on request

Coil Extension Tubing

Material: Polyurethane (PUR), reusable

20.0 cm (8"), max.
2.5 m (8 1/4 ft.) stretch length

Color	REF
red	20-10-442
blue	20-10-444
black	20-10-448

80.0 cm (32"), max.
10.0 m (33 ft.) stretch length

Color	REF
red	20-20-442
blue	20-20-444
black	20-20-448



Tubing


Material: Polyvinyl Chloride (PVC), for single use

7 mm O.D. / 3 mm I.D.
Length: per meter

Color	REF
red	20-20-652
blue	20-20-654

7 mm O.D. / 4.5 mm I.D.
Length: per meter

Color	REF
red	20-34-710-16
blue	20-30-710-16



Dual Lumen Tubing

Material: Polyvinyl Chloride (PVC), for single use

Length: per meter

Color		REF
red	6.0 mm O.D. / 3.0 mm I.D. and 6.0 mm O.D. / 3.0 mm I.D.	20-20-542
blue	6.0 mm O.D. / 3.0 mm I.D. and 6.0 mm O.D. / 3.0 mm I.D.	20-20-544

