

# Tourniquet Cuffs

## Product details

### Silicone Bladder

- Silicone material is elastic and acts as a low pressure cuff (as wider the cuff as lower the inflation pressures needed)
- Durable and high quality material
- Latex free

### Cuff ties

- Provide easy cuff positioning and security
- Color coded for identification of different sizes

### Contour cuffs

- Ideal for extremely conical shaped limbs

### Silicone Tubing

- Kink resistant
- Never cracks
- Latex free

### Positive locking connectors (PLC)

- For easy attachment and detachment
- For safe and leak-free inflation

## Single Cuffs



### Fabric Cuffs

- Washable at 60°C
- Bladder and fabric sleeve can be replaced



### Silicone Cuffs

- Autoclavable at 134°C (273°F)
- Easy to clean and dry quickly



### Silicone/Velcro Cuffs

- Autoclavable at 134°C (273°F)
- Velcro sleeve can be replaced
- Easy to clean and dry quickly



### Dispo Cuff

- Disposable, sterile peel pouch
- Low profile, welded cuff (no bladder)
- Eliminate risk of cross contamination



### Slim Cuff


- Washable at 60°C
- Low profile, welded cuff (no bladder)

Patient, Size	Fabric Cuffs	Silicone Cuffs	Silicone/Velcro Cuffs	Dispo Cuff	Slim Cuff
<b>Infant</b> 20 cm (8")	20-54-700	—	20-64-700	20-34-700SLZ-1	20-74-700
<b>Child</b> 30 cm (12")	20-54-710	20-94-710	20-64-710	20-34-710SLZ-1	20-74-710
<b>Arm</b> 35 cm (14")	20-54-711	20-94-711	20-64-711	20-34-711SLZ-1	20-74-711
<b>Arm,</b> 46 cm (18")	20-54-712	20-94-712	20-64-712	20-34-712SLZ-1	20-74-712
<b>lower Leg</b> 46 cm (18")	20-54-512 contour	—	—	20-34-715SLZ-1 contour	20-74-715 contour
<b>Leg</b> 61 cm (24")	straight 20-54-722 contour 20-54-522	20-94-722	20-64-722	20-34-722SLZ-1 contour	20-74-722 contour
<b>Leg</b> 76 cm (30")	straight 20-54-727 contour 20-54-527	20-94-727	20-64-727	20-34-727SLZ-1 contour	20-74-727 contour
<b>Leg</b> 86 cm (34")	straight 20-54-728 contour 20-54-528	20-94-728	20-64-728	20-34-728SLZ-1 contour	20-74-728 contour
<b>Leg</b> 107 cm (42")	20-54-729	20-94-729	20-64-729	20-34-729SLZ-1 contour	20-74-729 contour

## Plastic Cover

### Tourniquet Cuff Protector

The plastic covers keep the cuff free from soiling and allow sealing the cuff with adhesive tape without damaging the cuff.



REF	Size	Box
REF 20-11-111	small	50
REF 20-11-111-1	ditto but sterile	50
REF 20-11-333	medium	50
REF 20-11-333-1	ditto but sterile	50
REF 20-11-222	large	50
REF 20-11-222-1	ditto but sterile	50

## Measuring Tape

color coded, for easy identification of patient extremity and cuff size, reusable

Box 10      REF 20-95-111

ditto but for single use  
(Box of 1 block with 50 tapes)

Box 1      REF 20-95-222



## Double Cuffs



### Fabric Cuffs

- Washable at 60°C
- Bladder and fabric sleeve can be replaced



### Silicone Cuffs

- Autoclavable at 134°C (273°F)
- Easy to clean and dry quickly



### Silicone/Velcro Cuffs

- Autoclavable at 134°C (273°F)
- Velcro sleeve can be replaced
- Easy to clean and dry quickly



### Dispo Cuff

- Disposable, sterile peel pouch
- Low profile, welded cuff (no bladder)
- Eliminate risk of cross contamination

Patient, Size	Fabric Cuffs	Silicone Cuffs	Silicone/Velcro Cuffs	Dispo Cuff
<b>Infant</b> 20 cm (8")	20-50-700	—	20-60-700	—
<b>Child</b> 30 cm (12")	20-50-710	20-90-710	20-60-710	20-30-710SLZ-1
<b>Arm</b> 35 cm (14")	20-50-711	20-90-711	20-60-711	—
<b>Arm,</b> 46 cm (18")	20-50-712	20-90-712	20-60-712	20-30-712SLZ-1
<b>Leg</b> 61 cm (24")	20-50-722	20-90-722	20-60-722	20-30-722SLZ-1
<b>Leg</b> 76 cm (30")	20-50-727	—	20-60-727	—
<b>Leg</b> 86 cm (34")	20-50-728	—	—	—
<b>Leg</b> 107 cm (42")	20-50-729	—	—	—

