

Automatic Pressure Infusors

For irrigation e.g. during Arthroscopy, Laparoscopy and Hysteroscopy

Product Details

Precision Pressure Regulator

- Very precise pressure regulation
- maintains selected pressure and guarantees automatic compensation
- limits pressure to max. 300 mmHg

Protection Cover

- Protects the device from liquid and damage

Coil Tubing

- Built into Pressure Infusor
- Leak proof
- With male Luer Lock Connector

Universal Clamp

- Horizontal or vertical position for attachment to drip stand or rail

Pressure Infusor 100

with single channel

for compressed air

supplied with 4 m (13ft.) pressure hose

with universal clamp

REF 52-13-100

with stand and basket

REF 52-15-100



Pressure Infusor 200

with switch for the connection of two Pressure Infusion Cuffs

for compressed air

supplied with 4 m (13ft.) pressure hose

with universal clamp

REF 52-13-200

with stand and basket

REF 52-15-200



Pressure Infusor 200 ELC

with switch for the connection of two Pressure Infusion Cuffs

for electricity, available in 230, 110 or 100 VAC

with universal clamp

REF 52-13-200 ELC

with stand and basket

REF 52-15-200 ELC

