

| | | |
|---|---|--|
| VBM Medizintechnik Quality assurance | Summary of latex containing VBM-products | Date: 20.09.2011 Issue: 17 File: see below Page: 1 of 1 |
|---|---|--|

To keep our customers aware of any possible or potential reactions and/or side effects from the use of our products we provide the following list of devices and device components that contain natural latex..



Warning:
Following medical devices contain latex. Contact with these products may cause allergic reactions in some individuals.

All other VBM products are latex free according to current knowledge.

| Product | Article Nr. | latex containing component |
|---|--|-----------------------------------|
| Pressure Infusion Cuffs Spare Part Rubber Bladder | 50-xx-100 50-xx-500 51-xx-025 51-xx-050 51-91-055 52-xx-100 52-xx-150 52-xx-300 52-xx-330 52-61-335 52-xx-500 53-xx-100 53-xx-110 53-xx-150 53-xx-200 53-xx-300 53-xx-330 53-xx-500 | Rubber bladder including tubings |
| Guedel Airway with elastic Necktape (only size 0, 00, 000) | 40-26-007-1 40-26-008-1 40-26-009-1 40-06-040 | Elastic Tape |
| Endotape Oral/Nasal | 40-8018-044 | Elastic Tape |
| Tourniquets | 20-30-000 20-32-000 | Elastic Tape |
| Warming cushion ^{1.)} | WK-400x (xx) | Elastic Tape |

^{1.)} not classified as medical device