

Automatic Tourniquet 6000

Low-cost electronic Tourniquet without battery

- Single Channel
- Power: 100-240 VAC
- Lightweight Design: 1.7 kg (3.75 lb)

Automatic Timer

- Automatic time recording after cuff inflation
- Provides elapsed time after cuff deflation
- Audible signal every 30 minutes

Carrying Handle

- Easy to transport

Digital Display

- Large LED display is easy to read
- Precise monitoring shows actual cuff pressure

Regulation

- Precise inflation and deflation in increments of 5 mmHg
- Limits pressure to max. 600 mmHg



Flush Button

- To check for bleeding after surgery

Alarm System

- Visual and audible alarm system to detect leaks and pressure deviations

Coil Extension Hose

- Kink resistant
- Extendable, no hoses on the floor
- fixed to the machine, does not get lost

Positive Locking Connector

- Easy attachment and detachment
- For safe and leak-free inflation

Order information

Table Unit with
Carrying Handle
REF 14-12-600

Table Unit with Carrying
Handle on mobile stand
with Basket
REF 14-22-600



VBM Medizintechnik GmbH

Einsteinstrasse 1 | DE-72172 Sulz a.N. | Germany

Tel.: +49 7454 / 959610 | Fax: +49 7454 / 959633 | e-mail: info@vbm-medical.de | www.vbm-medical.de

