


# VBM Tourniquet Cuffs, conical

VBM offers a variety of contour cuffs. Contour cuffs have an arced design with a smaller diameter distally than proximally when applied. They enhance patient comfort, ensure a snug fit for conically shaped limbs, reduce the risk of mechanical shearing and guarantee an effective limb occlusion.

## Fabric Cuffs

- Washable and reusable
- Silicone Bladder and fabric sleeve can be replaced

Limb	Size	REF	Box
Arm/lower leg	46cm (18")	20-54-512	1
Leg	61cm (24")	20-54-522	1
Leg	76cm (30")	20-54-527	1
Leg	86cm (34")	20-54-528	1



## Disposable Cuffs

- Sterile peel pouch
- Low profile, welded cuff (no bladder)
- Eliminate risk of cross contamination

Limb	Size	REF	Box
Leg	61cm (24")	20-34-722SLZ-1	10
Leg	76cm (30")	20-34-727SLZ-1	10
Leg	86cm (34")	20-34-728SLZ-1	10
Leg	107cm (42")	20-34-729SLZ-1	10



**VBM** Medizintechnik GmbH

Einsteinstrasse 1 • DE-72172 Sulz a.N. • Germany

Tel.: +497454 / 959610 • Fax: +497454 / 959633 • e-mail: info@vbm-medical.de • www.vbm-medical.de

