

# VBM Slim Cuff

The new VBM Slim Cuffs complete the existing wide range of VBM Tourniquet Cuffs. The welded cuff offers a very low profile and is therefore lightweight and easy to apply.



- material: Nylon, easy to clean for multiple use, washable at 60°C (140°F)
- color-coded fixation tapes provide easy cuff positioning
- kink resistant silicone hose with positive locking connector
- inexpensive alternative to fabric- and silicone cuffs

## contour shape

- delivers uniform pressure
- fits limb of different size and shape
- roll resistant due to inside stiffener

Single Cuff straight shape		
REF	Application	Size
20-74-700	Baby	20 x 5 cm (8")
20-74-710	Child	30 x 6 cm (12")
20-74-711	Arm	35 x 8 cm (14")
20-74-712	Arm, long	46 x 10 cm (18")

Single Cuff contour shape		
REF	Application	Size
20-74-715	Arm/ Lower Leg	46 x 12 cm (18")
20-74-722	Leg	61 x 12 cm (24")
20-74-727	Leg, long	76 x 12 cm (30")
20-74-728	Leg, extra long	86 x 12 cm (34")
20-74-729	Leg, super long	107 x 12 cm (42")

**VBM** Medizintechnik GmbH

Einsteinstrasse 1 | DE-72172 Sulz a.N. | Germany  
 Tel.: +49 7454 / 959610 | Fax: +49 7454 / 959633 | e-mail: info@vbm-medical.de | www.vbm-medical.de

