

VBM MANUJET III

The Manujet III with Jet Ventilation Catheters is a necessity for emergency rooms, crash carts, ambulances and operating rooms as it guarantees a quick and efficient oxygenation of a patient

Emergency use

- As a life saving manoeuvre in the “cannot intubate – cannot ventilate“ situation for oxygenation to avoid a severe desaturation of the patient
- For prehospital use if there is an obstruction of the upper airway. Trans Tracheal Jet Ventilation (TTJV) is faster and simpler than a tracheostomy in emergencies with fewer complications

Elective use


- Micro laryngeal surgery
- Rigid bronchoscopy
- ENT surgery
- To assist a difficult fiberoptic intubation



MANUJET III

with 4m pressure hose, Jet Ventilation Catheters acc to Ravussin (13G, 14G, 16G), 100cm connecting tube, Endojet Adapter (for ET Tube) and Bronchoscope Adapter


REF
30-01-003



VBM Jet Ventilation Catheter

acc. to Ravussin

Allows a safe and atraumatic transtracheal access for oxygenation/ventilation with a manual jet (Manujet III) or an automatic High Frequency Jet Ventilator



Ordering information

ANESTHESIA SET

Jet Ventilation Catheter acc. to Ravussin
sterile packed, 5 pcs. per box

Size	Patient	REF
13G	Adult	30-02-013-1
14G	Child	30-02-014-1
16G	Infant	30-02-018-1

EMERGENCY SET

Jet Ventilation Catheter acc. to Ravussin
sterile packed, 1 piece per box with syringe

Size	Patient	REF
13G	Adult	30-02-913-1
14G	Child	30-02-914-1
16G	Infant	30-02-918-1

Connecting Tube, 100cm long with Luer Lock

REF
54-05-111

Endojet Adapter, for connection to the Endotracheal Tube, 15mm I.D. swivelling, 22mm O.D. with expiration holes, 50cm adjustable catheter with Luer Lock, sterile packed

30-03-111-1

Bronchoscope Adapter, for connection to the rigid bronchoscope, 15mm I.D. swivelling, with Luer Lock

60-70-005



VBM Medizintechnik GmbH

Einsteinstrasse 1 • DE-72172 Sulz a.N. • Germany

Tel.: +49 7454 / 959610 • Fax: +49 7454 / 959633 • e-mail: info@vbm-medical.de • www.vbm-medical.de

