

VBM Bi-Mask

Flexible Mask due to unique design with the combination of two materials

with silicone lip



silicone lip
 • anatomical shape, for ideal seal

handle
 • for a better grip
 • smooth surface for attachment of hook ring

silicone
 • transparent, clear view

plastic shell
 • molded into silicone for improved stability
 • color coded

with inflatable cushion



inflatable cushion
 • anatomical shape, for ideal seal

Bi-Mask with silicone lip

latex free, autoclavable up to 134°C

Size	Patient	Color	REF	Box
3	Adult small	yellow	35-65-223	1
4	Adult medium	red	35-65-224	1
5	Adult large	purple	35-65-225	1
6	Adult extra large	blue	35-65-226	1



Bi-Mask with inflatable cushion

latex free, autoclavable up to 134°C

Size	Patient	Color	REF	Box
3	Adult small	yellow	35-75-223	1
4	Adult medium	red	35-75-224	1
5/6	Adult large	purple	35-75-225	1



Hook Ring

Ø 30-33 mm,
made of Polysulfone (PSU)

REF 35-60-266
Box 10



Fixation Tape


latex free, for single use,
made of rubber

REF 35-30-255
Box 10

Harness

latex free, autoclavable
made of silicone

REF 35-70-255
Box 10



VBM Medizintechnik GmbH

Einsteinstrasse 1 | 72172 Sulz a.N. | Germany

Tel.: +49 7454 / 95 96 10 | Fax: +49 7454 / 95 96 33 | e-mail: info@vbm-medical.de | www.vbm-medical.de



CE
0123