


VBM Bi-Mask

Flexible mask due to the unique design with the combination of two materials. The silicone material of the mask is soft and comfortable for the patient. The plastic shell is molded into the silicone and gives the mask stability allowing the user to manipulate the mask for a perfect seal.



Bi-Mask

latex-free, autoclavable up to 134°C

Size	Patient	Color	REF	Box
3	Adult small	yellow	35-65-223	1
4	Adult medium	red	35-65-224	1
5	Adult large	purple	35-65-225	1
6	Adult extra large	blue	35-65-226	1



Hook Ring

Ø 30-33 mm, made of polysulfone (PSU)


REF 35-60-266
Box 10



Harness

latex free, autoclavable up to 134°C
Material: Silicone

Patient	REF	Box
Child	35-70-155	10
Adult	35-70-255	10



VBM Medizintechnik GmbH

Einsteinstrasse 1 • DE-72172 Sulz a.N. • Germany

Tel.: +49 7454 / 959610 • Fax: +49 7454 / 959633 • e-mail: info@vbm-medical.de • www.vbm-medical.de

