


VBM Stylets, Introducers and Tube Exchangers

Introducer

for difficult intubation

for single use, sterile

- lumen-design, Universal Connector included
- rigid pre-formed tip
- graduation marks



Size	Length	for ET Tube	REF	Box
14 CH	65 cm	≥ 5.0 mm	33-50-650-1	5

METTS (Muallem ET Tube Stylet)

for difficult intubation

single use, sterile

- core made of metal, malleable (maintains curvature)
- flexible pre-formed tip
- graduation marks



Size	Length	for ET Tube	REF	Box
8 CH	40 cm	≥ 3.5 mm	33-08-400-1	5
12 CH	65 cm	≥ 5.0 mm	33-12-650-1	5
14 CH	65 cm	≥ 6.0 mm	33-14-650-1	5

METTI (Muallem ET Tube Introducer)

for difficult intubation and tube exchange

for single use, sterile

- core made of plastic, semi-rigid (flexible)
- flexible pre-formed tip
- graduation marks



Size	Length	for ET Tube	REF	Box
12 CH	80 cm	≥ 5.0 mm	33-12-800-1	5
14 CH	80 cm	≥ 6.0 mm	33-14-800-1	5

Tube Exchanger

for extubation and exchange of Endotracheal Tubes

for single use, sterile

- lumen-design, Universal Connector included
- graduation marks



Size	Length	for ET Tube	REF	Box
11 CH	80 cm	≥ 4.0 mm	33-45-800-1	5
14 CH	80 cm	≥ 5.0 mm	33-46-800-1	5
19 CH	80 cm	≥ 7.0 mm	33-47-800-1	5



Universal Connector

15 mm / Luer Lock for oxygenation, fits all sizes with lumen-design
Material: ABS, for single use

REF 33-60-005 Box 10

Intubation Catheter

for fiberoptically assisted intubation

for single use, sterile

- lumen-design, Universal Connector included
- graduation marks



Size	Length	for ET Tube	for Fiberoptic Scope	REF	Box
19 CH	56 cm	≥ 7.0 mm	≤ 4.0 mm	33-60-560-1	5

VBM Medizintechnik GmbH

Einsteinstrasse 1 | 72172 Sulz a.N. | Germany

Tel.: +49 7454 / 95 96 10 | Fax: +49 7454 / 95 96 33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

