


# VBM Rescuefix

## Oral tube fixation for emergencies with integrated biteblock

In emergency situations a quick and safe tube fixation is vital. The Rescuefix allows a fixation without the use of adhesives as the integrated clamp fixes the tube safely and quick.



### RESCUEFIX

Patient	Color code	for ET Tube	REF	Box
Child	yellow	5,0 - 6,5mm	40-9005-111	10
Adult	white	7,0 - 10,0mm	40-9005-333	10

**VBM** Medizintechnik GmbH

Einsteinstrasse 1 • DE-72172 Sulz a.N. • Germany

Tel.: +49 7454 / 959610 • Fax: +49 7454 / 959633 • e-mail: info@vbm-medical.de • www.vbm-medical.de

