

# VBM Infusor 500 / 1000

Transparent Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes.

## ■ transparent

- allows clear view on the fluid bag and the fluid level

## ■ Wrap-Around Design

- fluid bag can be replaced quickly and easily
- consistent pressure on the complete fluid bag

## ■ PU-Material

- easy to clean
- reusable
- durable
- latex free



### Stop Cock

- allows quick and complete deflation of the cuff

### large Hand Inflator

- for quick inflation of the cuff

## Ordering information

Cuff can be supplied with the following three options for pressure control:

### • Shock Resistant Precision Manometer

- exact pressure control
- overpressure proved
- shock resistant and long-lasting

VBM Infusor 500 ml  
REF 57-08-050

VBM Infusor 1000 ml  
REF 57-08-100



### • Low-Cost Precision Manometer

- exact pressure control
- overpressure proved

VBM Infusor 500 ml  
REF 57-06-050

VBM Infusor 1000 ml  
REF 57-06-100



### • Pressure Indicator

VBM Infusor 500 ml  
REF 57-02-050

VBM Infusor 1000 ml  
REF 57-02-100



**VBM** Medizintechnik GmbH

Einsteinstrasse 1 | DE-72172 Sulz a.N. | Germany

Tel.: +49 74 54 / 95 96 10 | Fax: +49 74 54 / 95 96 33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

