

# Accessories

**Hook Ring**  
made of metal, Ø 33mm  
REF 80-11-266



**Hook Ring**  
made of PSU, Ø 30-33mm  
REF 35-60-266



**PEEP Valve**, autoclavable up to 134°C  
REF  
2.5-10cmH<sub>2</sub>O, 22mm O.D., green 80-13-001  
2.5-10cmH<sub>2</sub>O, 30mm I.D., green 80-13-004  
5-20cmH<sub>2</sub>O, 30mm I.D., blue 80-13-002  
5-20cmH<sub>2</sub>O, 22mm O.D., blue 80-13-005

**PEEP Adapter**  
for connection of PEEP Valve,  
autoclavable up to 134°C  
REF 80-13-003

**Harness**, latex free, autoclavable up to 134°C

Patient	Material	REF	Box
Child	Silicone	35-70-155	10
Adult	Silicone	35-70-255	10




**Reservoir Valve**  
autoclavable up to 134°C  
REF 80-12-030



**Reservoir Bag**   
for single use

Size	REF	Box
2600ml	80-12-010	5
600ml	80-12-020	5



**O<sub>2</sub>-Tubing** made of PVC   
Length: 2m, for single use  
REF 85-17-000



**Extension Tubing with Connector**  
made of Silicone, autoclavable up to 134°C,  
for extension between Resuscitation Bag and  
Non-Rebreathing Valve  
REF 80-22-040



**Hook** made of Silicone  
autoclavable up to 134°C,  
fits all Resuscitation Bags  
REF 80-10-999




**Rail Clamp**  
with holder for  
Resuscitation Bag  
REF 79-025-00



**Intake Valve for Demand Valve**, complete  
Material: Aluminium, autoclavable up to 134°C

Connection of  
- Demand Valve (photo 1)  
Via additional adapter connection of  
- Reservoir Valve with Reservoir Bag and  
oxygen Tubing (photo 2)

REF 80-16-900



# VBM Silicone Resuscitation Bag

latex free, autoclavable up to 134°C



**Adult 2000ml with palm holder**  
- with 40cmH<sub>2</sub>O Pressure Relief REF 82-10-100  
- without Pressure Relief 82-10-101

**Child 500ml**  
- with 40cmH<sub>2</sub>O Pressure Relief REF 80-10-200  
- without Pressure Relief 80-10-201  
- with manometer for pressure control 80-10-251



**Infant 250ml**  
- with 40cmH<sub>2</sub>O Pressure Relief REF 80-10-300  
- without Pressure Relief 80-10-301  
- with manometer for pressure control 80-10-351



# Face Masks made of Silicone, latex free, autoclavable up to 134°C

## Air Cushion Face Mask

Size	Patient	Color	without Hook Ring REF	with Hook Ring made of metal REF	with Hook Ring made of PSU REF	Box
2	Child	green	35-70-222	35-71-222	35-72-222	1
3	Adult small	yellow	35-70-223	35-71-223	35-72-223	1
4	Adult medium	red	35-70-224	35-71-224	35-72-224	1
5	Adult large	purple	35-70-225	35-71-225	35-72-225	1
6	Adult extra large	blue	35-70-226	35-71-226	35-72-226	1



## Mask for Children

### Rendell Baker Mask

Size	Patient	REF	Box
0	Neonate	35-60-000	10
1	Infant	35-60-001	10
2	Child small	35-60-002	10
3	Child	35-60-003	10



### Mask, round

Size	Patient	REF	Box
0	Neonate	80-11-000	10
1	Infant	80-11-001	10
2	Child	80-11-002	10



## Mask for Adults

### Mask with plastic dome

Size	Patient	REF	Box
3	Adult small	80-11-003	10
4	Adult medium	80-11-004	10
5	Adult large	80-11-005	10



### Full Silicone Mask

Size	Patient	REF	Box
3	Teenager / Adult small	35-60-223	1
4	Adult medium	35-60-224	1
5	Adult large	35-60-225	1



# Complete Sets

## Resuscitation Bag Set for Adults

- supplied in a rigid transparent carrying case with
- 2000ml Silicone Resuscitation Bag with palm holder
  - Air Cushion Face Mask size 5
  - Reservoir Valve
  - 2600ml Reservoir Bag
  - 2m O<sub>2</sub>-Tubing

REF 82-10-190

## Resuscitation Bag Set for Children

- supplied in a rigid transparent carrying case with
- 500ml Silicone Resuscitation Bag
  - Air Cushion Face Mask size 3
  - Reservoir Valve
  - 2600ml Reservoir Bag
  - 2m O<sub>2</sub>-Tubing

REF 82-10-290

## Resuscitation Bag Set for Infants

- supplied in a rigid transparent carrying case with
- 250ml Silicone Resuscitation Bag
  - Mask, round, size 1
  - Reservoir Valve
  - 600ml Reservoir Bag
  - 2m O<sub>2</sub>-Tubing

REF 80-10-390

- Case for Child and Infant Resuscitation Bags  
Case for Adult Resuscitation Bag

REF 80-10-080  
REF 80-10-090



# Respiration Cases

## Respiration Case for Adults

- complete with
- 2000ml Silicone Resuscitation Bag with palm holder
  - Air Cushion Face Mask size 3, 4, 5
  - Airways size 3, 4, 5

REF 30-92-100

## Respiration Case for Children

- complete with
- 500ml Silicone Resuscitation Bag
  - Mask, round, size 1
  - Air Cushion Face Mask size 2, 3
  - Airways size 0, 1, 2

REF 30-92-200

## Respiration Case for Infants

- complete with
- 250ml Silicone Resuscitation Bag
  - Masks, round, size 0, 1
  - Air Cushion Face Mask size 2
  - Airways size 00, 0, 1

REF 30-92-300



Wall Fixation  
for Respiration Case  
REF 30-91-001

